

ZC-A200

H type



380mm

Lowest chair position

760mm

Highest chair position

≥150kg

Carrying capacity
of dental chair



- The practical super-low seat position 380mm serves the elder, child and handicapped well.
- Patient gets good lumbar support from the stable, durable and wide chair.



LED operation light

LED operation light has three-axis rotation, manual switch and infrared sensor switch, which can control the on-off and adjust the intensity. The handle can be disassembled, cleaned and disinfected. Vertical distance of the light arm is 650mm.



No-sticker integrated assistant control

1. Prevent the liquid from entering the circuit board through the gap of the sticker, improve the reliability of the circuit. Side handle makes operation easier.
2. The strong suction and saliva ejector can be turned off with a delay of 10 seconds to ensure that the dirt in tubes is completely drained.
3. The strong suction and saliva ejector handles are made of a combination of aluminum alloy and silicone, light weight and good sealing, can be removed for cleaning and disinfection. The flow can be adjusted by the handle flow switch.

Large doctor table

More space for materials and tools. Transparent cover prevents the erosion of the solutions and facilitate cleaning.



Fashion box

The design is streamlined and beautiful. The side cover is magnet type, easy to open for maintenance.



Rotary ceramic spittoon

Facilitates spitting by patients and facilitates four-handed operation.



Enclosed power source

According to the requirements of European safety regulations, the high-voltage part is enclosed, moistureproof and splash-proof, prevent electromagnetic interference and provide stable power output.

Multifunctional pedal switch

The pedal switch can control the speed of the handpiece according to the strength of the step. It integrates dental chair control, blowing air, handpiece water, flushing of the cuspidor and gargling water, which is convenient for doctors to use.



Low-voltage heater

Heater is heated by a DC 24V external heating method, which avoids changing the water quality and heating effect due to the scale adhering to the heating rod and prevent dry burning.



Eco-friendly PU Cushion | Standard



Continental Leather Cushion | Option

